
MILKYWAY PULSAR LIGHT PAD

NUOVISSIMA TECNOLOGIA LED A 3 CHIP

USER MANUAL

Before starting your first session, please make sure to read this manual as it explains important operating procedures.

Notes: Each LED module has three chips: 1pc **670** nm, 1 pc **830** nm and 1 pc **850** nm.

The light mat uses a combination of red 670nm and infrared 830/850nm LEDs, which are chosen to give the ideal wavelengths with a biological activity found in the scientific research of the last three decades.

Specifications:

Power: **170W**

Size: **200*70 cm**

Led quantity: **1620** pcs

G.W: **4.5 Kg**

Wavelengths: **670-830-850 nm=1:1:1**

- **670nm Red Light**

This light can be seen by the human eyes and targets skin cells.

- **830-850nm Near-Infrared light**

This light can't be seen by human eyes and targets deep tissues.

Benefits backed by science:

Research shows that these wavelengths trigger an optimal biological response to:

- Eliminate spots and wrinkles
- Stimulate collagen synthesis (smooth and firm skin)
- Improve skin tone and texture
- Relieve muscle and joint pain
- Improve blood circulation
- Help with weight loss
- Speed up wound healing
- Exert anti-inflammatory action
- Help for hair regrowth

- **FIRST MODE (ALL – 670/830/850 nm)**
- **SECOND MODE: ONLY RED (670 nm)**
- **THIRD MODE: 3-5 -> 830 nm AND 850 nm TOGETHER**

PLEASE NOTE INFRARED LIGHT IS NOT VISIBLE TO HUMAN EYE SO THE RELAVITE LED CHIPS SEEM “OFF”.

- Photobiomodulation can induce significant improvement in skin complexion, tone, elasticity, roughness, softness, smoothness, firmness and collagen density. Studies have shown efficacy and safety for skin rejuvenation and increase in intradermal collagen and elastic fibers, with activation of fibroblasts.
- Photobiomodulation is successfully used for healing wounds, burns, amputations, skin grafts and other injuries.
- Photobiomodulation stimulates the production of energy (ATP) and is therefore important for giving vigor and vitality and increasing physical and mental performance.
- Photobiomodulation improves cell survival, stimulates stem cells and protein synthesis.

- Photobiomodulation increases antioxidant defenses and reduces oxidative stress. One of the most observed effects is a reduction of inflammation.
- Red/infrared light can improve circulation due to its vasodilating abilities.
- Photobiomodulation stimulates the synthesis of melatonin, one the most amazing molecules in the body, which is a powerful antioxidant and regulates circadian rhythms.

SUGGESTED USE:

1. Ideal treatments last 30 min, but the exposure time can be adjusted according to personal features and needs: for people that have inflammatory conditions and/or need energy boosts, it is possible to increase the exposure. On the contrary, energetic and hyperactive people may need less time, and generally no more than the suggested average of 30 min.
2. For superficial skincare it is possible to use the red light only or in combination with near infrared; for longer treatments, for sensitive people and late in the afternoon (especially for people with insomnia) it is possible to use a lower power output and/or near infrared light only.
3. For deep skincare and for nourishing organs, bones muscles, spinal cord, etc. it is suggested using the mode "ALL" and "P5".
4. Using the pulse mode (10 and 40 Hz) allows to penetrate even deeper to nourish tissues and sooth inflammation; it is also really good for elderly people.
5. We suggest alternating continuous and pulsed light mode between different treatments or within the same treatment in order to provide a wider range of stimuli.
6. For people affected by seizures or any condition that could make prone to be hypersensitive to pulsed light, it is better to use continuous light mode or switching off red light and use near infrared light only.
7. "H" mode is suitable for improving mood and facilitate relaxation, so it can be used with a lower power.

In general, it is best to avoid using pulsed light and maximum power output in the late afternoon.

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